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*Mustoe TA et al. Plast Reconstructive Surg 2002;

110:560-71

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Kelo-cote® Spray is a lightweight, self-drying silicone gel spray for the treatment of scars. Kelo-cote® Spray rapidly dries to form a sheet; this Silicone Gel Sheet layer is gas permeable, flexible and waterproof.

Kelo-cote® Spray binds to the stratum corneum (the outer layer of dead skin cells) forming a protective barrier against chemical, physical and microbial invasion of the scar site while assisting with hydration. This creates an environment which allows the scar to mature through reduced collagen synthesis cycles, and improves the physiological and cosmetic appearance of the scar.

Silicone gels have been shown to be effective in treating abnormal scarring in the form of hypertrophic scars and keloids. The 2002 International Clinical Recommendations on Scar Management recommended “silicone gel should be the first line of treatment in the initial management of scars and in the prevention of hypertrophic scars and keloids.”*

Kelo-cote® Spray can be used on joints, flexures, hairy areas and large area scars such as those caused from burns or major trauma.

Intended use

Kelo-cote® Spray is intended for the management of old and new scars, including hypertrophic and keloid scars, resulting from general surgical procedures, trauma, wounds and burns.

Directions for use

- Shake the can well before use.
- Ensure that the affected area is clean and dry.
- Remove cap from the can.
- Hold the can 15-20 cm from surface to be treated.

- Spray along the length of the scar in a light sweeping motion.
- Kelo-cote® Spray should be applied in a very thin, even coating with no holes or gaps on the affected area.
- Allow to air dry before covering with sunscreen or pressure garments.
- Kelo-cote® Spray should dry in 2-3 minutes. If the Kelo-cote® Spray is pooling, dripping or running you have sprayed too much. Gently remove the excess and allow the drying process to continue.

Kelo-cote® Spray should be applied twice daily for 90 days for maximum effect.

Can be applied under pressure garments and sunscreen.

Recommended minimum treatment time is 60–90 days.

It may take longer than 90 days for older and larger scars to achieve full benefit. Extended use is recommended if continued improvement is seen.

Caution

- Kelo-cote® Spray contains flammable gases to help deliver the silicone.
- Extremely Flammable. Keep away from heat.
- Do not spray on open flame or incandescent material.
- Keep away from sources of ignition.
- Do not smoke while using the product.
- Pressurized container: Protect from sunlight and do not expose to temperatures exceeding 50°C.
- Do not pierce or burn can, even after use.
- Do not inhale Kelo-cote® Spray.

Do not use Kelo-cote® Spray near the eyes, nose or mouth.

Avoid direct contact with mucous membranes, third

degree burns and open wounds.

For external use only. If irritation occurs, discontinue use and consult your physician. Keep this product out of the reach of children.

Store below 25°C.

Do not use after the expiration (EXP) date printed on the can .

The lot number is printed on the bottom of the can.

Ingredients

Butane, Pentane, Polysiloxanes, Silicon Dioxide

Medical device registration number: 2140001

Manufacturer:

Advanced Bio-Technologies, Inc.

Suwanee, GA 30024, USA

www.kelocote.com

US Patent 5,741,509

International Patents

Registration holder:

Rafa Laboratories Ltd., P.O.Box 405, Jerusalem

9100301 Tel. 02-5893939

*Mustoe TA et al. Plast Reconstructive Surg 2002;

110:560-71

February 2012

Kelo-cote Spray PIL PB0214-05